


Specification of RF AF297MO-2019- -KBL

01: Negatoscope (X-ray wall viewer):



Body: 0, 70 mm sheet

Diffuser: Opal PMMA

Ballast: EN 60921-EN 61347-2-8-EN 61347-2-3

Starter: EN 60155

Base: EN 60400

Cable: 1 mm NYA

Coating: Electrostatic powder coating

General information

X-ray wall-mounted viewer or illuminator; used to visualize developed X-ray films in an X-ray department or operating theatre.

Specifications

Specifications:	Designed to read X-rays; it is classified as a radiology accessory.
	3 essential characteristics must be checked to ensure the quality of a Negatoscope X-ray wall-mounted film viewer:

	1) Coloured temperature tubes: minimum 3600 to 4000 Kelvin , "daylight" fluorescent tubes.
	2) Homogeneity of the illuminated surface: ensure that all zones are adequately illuminated.
	The maximum tolerated deviation between the centre and the corners is a variant of 30% (i.e., 8cm from the edges of the illuminated surface).
	3) Light power: The minimum luminance for a Negatoscope with an on/off switch must be set greater than 1700cd/m ² .
	Manual processing: a drip tray is recommended on the light box..
	Size varies slightly dependent on manufacturer.

02- BP set (Sphygmomanometer and stethoscope) Stand type



Technical specification:

Handwritten signature

Handwritten signature


Two-piece stand column can be adjusted 80 to 125cm
mounted with 5pcs castors.

Standard adult 2 tube nylon cuff

2m spring tube

Measure scope: 0-300mmHg

03- Pulse oximeter (infant):




Technical specification:

- The most dependable oxygen saturation sensor for accurate readings.
- It combines the convenience and performance of a disposable sensor with the benefits of a reusable sensor. The reusable flex sensor is made of durable silicone and secured with a durable flexi wrap adhesive that is also disposable. This helps you keep a check on your child's health without disturbing their routine. This pulse oximeter probe for infants is great for extended monitoring and sleep studies.
- Packaged with 25 flexi wraps
- Accurate and dependable
- SpO2 readings start with Pure Light Sensors
- 1 Meter, 2 to 20 Kg (4.4 to 44 lbs)

[Handwritten signatures and initials in blue ink]

04-Nebulizer machine:



Specification

Portable, compressor driven nebulizer pump

Easy to operate and sturdy design for use in demanding environment

Cup capacity: m

in 3 ml to max 10 ml

Nebulization, approx: 0.3 to 0.8 ml/min

Residual volume: max 2 ml

Flow delivery range (under load): 0.5 to max 10 l/min

Mass median aerodynamic diameter of particles (MMAD): 0.5 to 4 μ m

Operating temperature range: +5 to +50 C

Operating relative humidity range: up to 65 %

Compressor requires no maintenance and delivers oil-free flow


Device, with carry handle, is provided with accessory storage compartment

Power supply: 220 V / 50 Hz

Power consumption, approx: 80 W

[Handwritten signature]
[Handwritten signature]

05-Phototraphy Machine blue light



[Handwritten signature]

[Handwritten date: Jan 19, 1961]


Parameter	Specification
Power	110 - 230Vac, 50/60Hz, 100VA
Dimensions	588 x 244 x 110mm (WxDxH) 5kg (excluding stand)
Light Output	Wavelength: 440 – 470nm
Regular Mode	25 μ W/cm ² .nm (mean)
Intensive Mode	>35 μ W/cm ² .nm (mean)
Therapeutic Mode	10 – 100 μ W/cm ² .nm, steps of 5 μ W/cm ² .nm (peak)
Illuminated Area	600 x 300mm at 360mm distance
LED life	>60,000 hours
Temperature Monitor	
Range	25 - 43°C
Accuracy	±0.3 °C
Alarms	High & Low Temperature, user programmable
Mobile Stand	Optional Accessory
Height	1200 – 1550mm, adjustable
Base	750 x 707 x 100mm (WxDxH)
Weight	16kg
Compliance	European Medical Device Directive 93/42/EEC EN 60601-1; EN 60601-1-2; EN 60601-2-50
Compliance	EN1400; US 16-CFR-1511

Baby receiving blue light phototherapy, or bili therapy, to treat neonatal jaundice (hyperbilirubinemia). Jaundice is a condition in which there is an excess of the yellow pigment bilirubin in the body's tissues, seen as a yellow discoloration of the skin and eyes. Bilirubin is the breakdown product of haem, the oxygen-carrying chemical in red blood cells. In neonatal jaundice, failure of the liver to process the bilirubin leads to its accumulation in tissues. This is very common in newborn babies, and is usually harmless. The blue light breaks down the pigment, allowing it to be excreted naturally. The baby's eyes are covered to prevent any damage to the retina.

[Handwritten signature]

[Handwritten signature]

06- Resuscitation table:



TECHNICAL SPECIFICATIONS:

Mobile resuscitation table for neonatal patients, with warmer

Sturdy construction, with 4 anti-static swivel castors and 2 brakes

Side handle for easy moving

Integrated support for 2 x 10 l oxygen cylinders

Side rail and IV pole for mounting accessories and equipment

With lower drawer for storage

Side panels, transparent acrylic resin, collapsible/foldable for complete accessibility

Vertical column with controls and display


[Handwritten signature]

[Handwritten signature]

Table:

surface 60 W x 76 L cm (minimum guaranteed dimensions)
mattress covering the entire table surface, thickness 5 cm
mattress cover flame-retardant, waterproof and cleanable with hospital-grade disinfection products
Heating system:
overhead radiant heater, fixed height
radiant heating by quartz/ceramic elements, parabolic reflector and safety grid

07- Complete blood count machine (CBC hematology Analyzer):



Technical specification:


- Use a blood sample to classify and quantify 12 different blood characteristics
- It contains safety lockouts, fail-safe and alert mechanisms
- Principle: Impedance for WBC, RBC, PLT counting
- Throughput: 60 samples per hour
- 3-part differentiation of WBC, 20 parameters +3 histograms
- 3 counting modes: whole blood, capillary blood and prediluted
- Micro sample volume: 10 μ L
- Storage: up to 100,000 sample results (including histograms)

Handwritten signature

Handwritten signature

- Support external printer, barcode and LIS

09: Microscope:



Technical specification:

Optics	
Lighting	Transmitted light
Lamp type	LED
Magnification	40x, 100x 400x, 1000x
Lens type	Plan, achromatic
Image scale	4, 10, 40, 100
Condensor	Abbe (Securely mounted)

Handwritten signature and date: 15/11/2017

Power supply	Power plug
Optical system	infinity
Objective 1	Plan 4x, NA = 0.1, AA = 27.8 mm, DK = -
Objective 2	Plan 10x, NA = 0.25, AA = 8 mm, DK = -
Objective 3	Plan 40x, NA = 0.65, AA = 0.60mm, DK = 0.17
Objective 4	Plan 100x, NA = 1.25, AA = 0.13 mm, DK = -
Capally brightfield	yes
Mechanics	
Type of build	binocular
Viewing posture (°)	30° angled eyepiece
Focusing	Coarse and fine movement
Revolver	4
Special features	
Built-in digital camera	-
Rechargeable battery operation	-
Filter holder	yes
Luminous field diaphragm	yes
Dust cover bag	yes
Centering stage	yes
General	
Length (mm)	238
Width (mm)	197
Height (mm)	385



Weight (g) 6000


Series CX23

Colour white

Field of application

Biology and Medicine yes


10-Plaster cutter large size:




Handwritten signature and scribbles in blue ink.

TENS Electrodes - Value Wired 2x4 Replacement Pads for TENS Units - 20 TENS Unit Electrodes -
2"x4" Wired TENS Unit Pads- Discount

11-plaster cutter:



12-baby warmer:



Handwritten signature and date:
11/11/11

Technical specification:

(Indicative specifications)

Standalone system with height adjustment and can be conveniently placed over any baby cradle or trolley.

Micro Processor Controlled with Servo & Manual modes of operation.

This infant warmer provides automatic over heat and over temperature protection in both manual and servo modes.

Quartz Infrared heater element with very effective heat Radiation

Thermistor based high sensitive Temperature Probe ensures an accuracy of $\pm 0.2^{\circ}\text{C}$. Probe is inter-changeable without calibration.

The Parabolic Reflector is designed to send uniform heat to the bed surface.

Memory Backup to restore set temperature and previous mode, in case of power failure.

Swiveling Capacity of the Heater Unit for taking X rays.

Displays set temperature & infant temperature simultaneously.

Comprehensive Alarms and safety cut-offs.


Wide base for stability and Anti Castor wheels with brakes for good mobility.

Easy to use and clean.

IV Pole, Infusion Pump Mounting Pole, Tray for Vital Sign Equipments, Halogen Lamp for observations.

Delivered with Baby Trolley for use with free standing warmer.

135-Steel catheter different size:




Handwritten signatures and initials in blue ink.

Item specifics:

- Chastity Device Type: Catheters & Sounds
- Model Number: YYGM
- Size: 9 Different PCS
- Material: Stainless Steel
- Brand Name: IWANT
- Stainless Steel. 9 Different Pcs. Secret delivery, nobody knows what's inside.
- Quantity: 1 Lot (9 Piece / Lot)
- Package Size: 30.0 * 10.0 * 6.0 (cm)
- Gross Weight/Package: 0.9 (kg)

12-Annual speculum 3 sizes:



Technical specification:


Mathieu Rectal Speculum 90mm Length x 15mm 195mm PH675320

90mm in length and 15mm in diameter, our Mathieu Rectal Speculum is 195mm in length and is made from the highest medical grade materials

[Handwritten signature]

[Handwritten signature]

13-Ophthalmoscope:



Technical specification:

Rechargeable ophthalmoscope

1.CE and ISO certification

2.high quality, factory price, professional service

3.aluminum box

Material ABS Plastic


Usage Hospital

Power 1.75 W

Specification 300x200x100 mm

Weight 1 Kg

14-Autoclave large size Horizontal 90lit:



General information

Autoclave for hospital structures performing surgical activities, in outlying hospital or emergency surgical set-ups where running water and sources of energy are scarce.

Specifications

Materials:	Lid, chamber, bottom grid: stainless steel
	Gasket: silicone.
Specifications:	An autoclave is a strong container used for high-pressure steam sterilization.
	Manually-operated 90-litre single-chamber vertical or horizontal autoclave, designed to function on 3 different types of energy sources (gas, kerosene and electricity) and no running water. The energy sources can be used separately or in combination, according to the power supply situation
	The autoclave is delivered with a pressure kerosene burner.
	- Power supply:
	Electricity 220-380V
	Electrical system has a manual switch enabling it to be adapted to a 4500W, 3000W or 1500W power source.
	Kerosene or gas
	- Container:
	Volume: 70 liters dabble chamber
	Internal chamber size: 400x700mm

[Handwritten signatures and initials in blue ink]

	Valves: different color-coded gates, to facilitate procedures
	Steam exchanger: regulating valve with a counterweight of a 1.2 bar or 2.2 bar
	Drying system: surface condensation (condenses exhaust steam through contact with cold metallic surfaces - oil)
	It can run an autoclave without running water.
	Lid: manual locking system.
	Delivered with service and operators manual: English and French.
Packaging:	Unit presentation: Each autoclave is delivered with the following items:
	- 2 stainless steel sterilization baskets
	- 1 GASKET, for cover type OPA, since 2001
	- 1 HEATING ELEMENT, 4.5kW, 230/400V
	- 1 BURNER, KEROSENE, with foot pump
	- 2 DRUM, STERILIZING, 16cm, diameter 29cm, lateral ellipse
	- 2 DRUM, STERILIZING, 24cm, diameter 34cm, lateral ellipse
	- User and maintenance manual: French and English
	Serial number should be written on the outside packaging.
	Packed in a solid wooden case, waterproof.
	Weight of autoclave, packed: 240kg
	Size of autoclave, packed: 0.92x0.92x1.80m: 1.53m ³
Other:	Always mention the autoclave's serial number when ordering spare parts.


15- Ultrasound machine probe (convex):

1 x Convex abdominal probe, frequency range: 2.5 / 3.5 / 5.0 MHz

1 x Convex transvaginal probe, frequency range: 4.5 / 6.5 / 8.0 MHz




16-Blood Shaker machine:



Scientific Medical LED Digital 3D Blood Shaker Machine:

Specifications	SK-D1807-S
Tilt angle	7°
Max. load capacity(with platform)	3kg
Motor type	DC motor
Motor input power	16W
Motor output power	10W
Speed range	10~80rpm
Speed display	LED
Timer display	LED
Time settings range	1min~19h59min
Voltage	100~240V,50/60Hz
Power	20W
Weight	2.7kg
Permissible ambient temperature and humidity	5~40°C, 80%
Dimension[W×D×H]	270×330×130mm
Protection class	IP21

17-Hydrcolator Machine with 4 heating pack (lumber and cervical):




Technical specification:

- Extra-large tank and oversized drain.
- A second thermostat guards against low water level and overheating.
- Overall Dimensions 32"W x 20"D x 33"H.
- High quality stainless steel
- Thermostatically controlled temperature
- 3" (8 cm) swivel-type rubber casters for silent, friction-free movement of mobile units
- No plumbing required
- The heater is thematically controller and maintain in the unit at a temperature between 760C and 800C. Available Models: 4 Steams Pack & 8 Steam Packs.

18-Dental unite complete set:

[Handwritten signature]
[Handwritten signature]



Technical specification:


- 24V DC motor for chair(safety system) x1 set
- Auto spittoon flushing and cup filler control system x1 set
- 3 programs inter-lock control system x1 set
- Instrument tray with air break (down-mounted) x1 set
- Built-in floor box with main switch x1 set
- Rotatable hand piece holder x1 set
- Assistant system x1 set
- Glass spittoon (Removable) x1 set
- Double articulate headrest x1 set
- Round foot switch x1 set
- Seamless cushion(PU) ,4-way joystick x1 set
- High suction and saliva ejector system x1 set
- outer water tank x1 set
- LED-2 operating light with sensor x1 set
- LED film viewer x1 set
- 3-way syringe(hot/cold) x1 set
- Built-in floor box x1 set
- ST-301 doctor stool x1 set
- **Optional Accessories**
- External junction box x1 set
- ST-603 doctor stool x1 set
- Soft leather cushion x1 set
- Automatic cup filler with sensor x1 set

[Handwritten signature]

[Handwritten signature]

- TPC Multi-media endoscope image system x1 set
- Imported/domestic 4 hole high speed hand piece x1 set
- Imported/domestic 4 hole low speed hand piece x1 set
- TPC LED curing light x1 set
- ART magnetostriction scaler(Built-in) x1 set
- Assistant element, designed for either 2 hand or 4 hand treatments, is easy to reposition for unrestricted patient access. Equipped with up to 4 instruments, is always ready for all types of treatments. Access to patient is easy thanks to the double articulated assistant arm. Turning force is adjustable at two joints to optimize the movement for precise working positioning. Safety anti-collision system protects your chair operation during patient positioning.

19-Dental set



Handwritten signature
1
1, 2, 3, 4, 5

DATA SPECIFICATION

PRODUCT NAME	: DENTAL INSTRUMENT SET
BRAND / ORIGIN	: DIMEDA / GERMANY
TYPE / CATALOGUE No.	: VARIOUS

Consisting of :

Dental Set

DENTAL diagnostic set

1	48.367.00	ADAMS napkin holder	1 pc
2	82.702.22	mouth mirror plane size 4=22mm	24 pcs
3	82.728.01	mouth mirror handle	24 pcs
4	83.031.03	London College tweezer fig. 3	24 pcs
5	83.346.10	WILLIAM FLY meas pmha & expl	6 prs
6	83.391.17	explorer 1end/2p. fig. 17	24 pcs
7	83.391.09	explorer 1end/2p. fig. 9	24 pcs
8	83.391.23	explorer 1end/2p. fig. 23	24 pcs

DENTAL surgical set

9	06.103.00	Scalpel handle no. 3 12cm/4 1/4"	2 pcs
10	06.110.00	Scalpel blade sterile fig. 10	2 pcs
11	06.115.00	Scalpel blade sterile fig. 15	2 pcs
12	08.340.11	Iris scissors 11,5 cm	2 pcs
13	08.341.11	Iris scissors 11,5 cm	2 pcs
14	10.260.12	SEMKEN forceps serrated	2 pcs
15	10.270.12	SEMKEN Forceps, 1x2 teeth, 12,5cm	2 pcs
16	18.312.01	LANGENBECK retr. 21cm 30x11 mm	1 pc
17	18.312.04	LANGENBECK retr. 21cm 40x11 mm	1 pc
18	24.374.17	MATHIEU needleholder TC	2 pcs
19	32.101.14	FRIEDMAN rongeur 14 cm	1 pc
20	50.151.19	LUNIATSCHEK gauze packer	2 pcs
21	82.652.00	water syringe self fill, metal	1 pc
22	82.873.11	HEMINGWAY curette double f.1/1	2 pcs
23	82.873.22	HEMINGWAY curette double f.2/2	2 pcs

Dental Surgical Treatment


24	06.103.00	Scalpel handle no. 3 12cm/4 1/4"	1 pc
25	06.104.00	Scalpel handle no. 4	1 pc
26	06.110.00	Scalpel blade sterile fig. 10	1 pc
27	06.115.00	Scalpel blade sterile fig. 15	1 pc
28	06.120.00	Scalpel blade sterile fig. 20	1 pc
29	06.402.01	Scalpel metal handle fig. 1	1 pc
30	06.408.01	PARKER scalpel delicate fig. 1	1 pc
31	06.408.02	PARKER scalpel delicate fig. 2	2 pcs
32	08.130.14	DEAVER scissors str.bl/bl 14cm	2 pcs
33	08.824.12	WAGNER scissors	1 pc
34	08.825.12	WAGNER scissors	1 pc
35	08.830.19	TOMATI gum punch	1 pc
36	10.102.14	Forceps dressing 14,5 cm	2 pcs

Exclusive Agent : PT. EJA Sarya Perdana
 Graha Mas Pemuda Blok AD-12 Jl. Pemuda, Rawamangun Jakarta 13220
 Phone : +62 21 471 2269 fax : +62 21 471 2213
 e-mail : eja-sarya@centrin.net.id

Page : 1 of 4

[Handwritten signatures and initials in blue ink]

20-Fetal Doppler (TCG Machine):



Description

Fetal Doppler Heaco L6S allows you to monitor and diagnose the heart rate of the fetus in pregnant women. It works on the principle of ultrasound, but uses very little power, so you can use it at home, without harming either the mother or the child.

Fetal Doppler consists of a monitor with a built-in speaker and sensor. The sensor is connected to the unit via a cord and is mounted on the housing. Heaco L6S very small and portable, has dimensions 32h138h96 mm and weighs only 270 g heart rate is measured in the range of 50 - 210 beats / min, and is displayed on the LCD screen. The device operates on two AA batteries.

The device can be applied to the 16 th week of pregnancy and more. Calculating heart rate is performed by digital signal processing (DSP), which guarantees the accuracy and safety.

1. Fetal Doppler;
- 2 x AA batteries - 2 pcs .;
3. Branded bag;
4. Instructions;

5. Warranty;

6. A copy of the registration certificate of the Ministry of Health of Ukraine;

7. A copy of the Certificate of Compliance UkrSEPRO. Attributes fetal doppler HR Performance

Monitoring from 16 weeks

LED Display

The operating frequency of 1.0 MHz

Operating mode Continuous / pulsed Doppler

The strength of the ultrasound $\leq 5 \text{ mW / cm}^2$

Measuring range 50-210 beats / min

The accuracy of $\leq \pm 2 \text{ beats / min}$

Warranty 12 months

21- Coagulation machine (electrocuted Machine)



Technical Specification:

[Handwritten signature and initials in blue ink]

COAGULATION

Mixing by magnetic motor / stirrer INR / ISI Value % / Seconds Auto / Manual start

Parameters PT / APTT / Fibrinogen

Filter 405 nm / 630 nm.

ELISA STRIP READER

Wells 8 position

Filters 405 nm/450 nm / 630 nm

Methods Cut-off / Multi standard with EP / KIN / Fix Time

O.D. Range 0 - 3.5A

OPTICAL SYSTEM

Principle Absorption spectrometry

Range of measurement From - 0.1 to 3.0A

Nominal spectral range Automatic via filter wheel

Filters 340, 405, 450, 505, 546, 578, 630 and 670nm interference type

Bandwidth 10 nm (± 2 nm)

$\pm 5\%$ at 0.10 A

Reading error $\pm 2\%$ at 1.00 A

$\pm 2\%$ at 2.00 A

Light source Halogen lamp

Zero Automatic

Settings Monochromatic and Bichromatic readings.

Detector Schottkey type Photodiode - Silicon

THERMOSTAT SYSTEM

Peltier regulated 25°, 30° & 37°C programmed

Accuracy $\pm 0.2^\circ\text{C}$

CUVETTES - DUAL MODE

Flow cell of 18 μL

Removable cuvettes Macro, Semi-micro & Micro

PRINTER

Type Fast Thermal Graphic

Printing speed 50 mm/sec

Dots per line 384 dots/line

Paper Thermosensitive - 57 mm wide, Easy Loading

Programming On/Off test wise

SCREEN AND KEYBOARD

Screen Graphic LCD 240 x 128 pixels

Back Light Blue LED Back lit

Keyboard Tactile Soft - Touch

SUCTION SYSTEM

Continuous flow system with peristaltic gear pump.
Sipping volume can be programmed from 200 μ l to 5.0 ml.

PROGRAMMING

Table of tests	up to 235
Table of units	up to 9 units
Quality control storage of the last 30 results	
Memory	512 KB EEPROM


QUALITY CONTROL

2 Controls per test - Normal & Abnormal
Levy-Jennings charts
Westgard's rules

POWER / DIMENSIONS

Voltage	110-220 V
Frequency	50 - 60 Hz
Maximum power	110W
Temperature	10° to 40° C
Maximum relative humidity	80%
Dimensions	45 cm x 41 cm X 17
Weight	8 kgs

22- Micro hematocrit reader:




For determination of packed cell volume, using 75 mm capillary tubes with total fluid column heights between 40 and 70 mm

Blue
Handwritten notes and diagrams in blue ink, including a series of vertical lines and a small sketch of a tube or column.

23-Ear thermometer:


- Technical Specification:
- Accessories: 21 pieces of disposable hygienic probe covers, Protective cap, 1x 3V lithium battery, Child's safe battery door. ** Source: Ipsos, online study among n=802 GPs and Pediatricians'



- Ear thermometer.
- Accessories: 21 pieces of disposable hygienic probe covers, Protective cap, 1x 3V lithium battery, Child's safe battery door. ** Source: Ipsos, online study among n=802 GPs and Pediatricians in four European markets (n=200 each in GE, UK, FR, NL), fielded in Jan-Feb 2017.
- Fast reading.
- Recalls last reading.

24-Vision measuring Machine (Auto refractor):

Handwritten signature and scribbles in blue ink.




The Vision Measurement System is a portable vision measurement system offered as a cost effective solution for manual, small part inspection. Designed with the practicability in mind, the DM-WVL200 has a compact design, lens power supply by USB or external power, fine-tuning device, and an external computer. These features provide for high portability, convenience, and easy upgrades. Add OVM Software to handle measuring applications with ease.

Features:

- Pretty design with compacted frame, steady performance and facilitated operation
- Aluminum alloy design is light, better rigidity and high heat conductivity
- New Z-axis transmission device with precision fine-tuning and rapid move
- Step zoom lens doesn't need any linear calibration with software
- High resolution 130,000 pixels CMOS without image card
- USB2.0 interface supports plug and play, can be linked to notebook PC
- The power supply is by the PC USB interface or the changeover device

25-Head light:

Handwritten signatures in blue ink, likely for approval or verification.



Light source: white light

Color temperature: 6000K

Power: 1W

Voltage: 5V

Luminous flux: 180lm

Use one 18650 rechargeable lithium battery (not included)

Two buttons, independent switch

COB light source, three-speed adjustment

Four mode adjustments

USB charging, overcharge protection, safe to use

26-Folding screen:



Handwritten signature and date:
 11/11/15

Technical Specification:

Mobile Folding Screens are high-quality medical screens suitable for use in sports club first aid rooms, physio rooms, treatment areas, gyms, spas, hotels etc... Perfect for providing privacy to clients and patients. The screens are articulated so they fold down to a manageable size for storage and are easy to move around on their built-in wheels. Because this design features multiple leg supports, stability of the screen is not a problem.

- Panels are replaceable and manufactured from flame-retardant, wipe-clean PVC
- Swivelling castors allow free and easy movement
- Minimal assembly is required

Sunflower Mobile Folding Screens are available in 9 great designs, including plain colors and graphics. The Mr. Men screen will surely be a favorite in Rugby club treatment rooms. These screens are made to order in Great Britain.

Three panel screen measures: 172cm(H) x 188cm(W) x 3cm(D) x 38cm (across legs).

Four panel screen measures: 172cm(H) x 250cm(W) x 3cm(D) x 38cm (across legs).

Five panel screen measures: 172cm(H) x 312cm(W) x 3cm(D) x 38cm (across legs)


27- Centrifuge electrical:



Weight (kg)	8
Warranty	one year
Motor Type	Universal Motor
Brand	Sigma
Portable	Portable

[Handwritten signatures and initials in blue ink]

28- SMR set (DNS)




TEM DETAILS:

- 1 B.P Handle No.3
- 1 Freer Dissector & Elevator
- 1 Freer Knife Str.
- 1 Freer Chisel
- 1 Killin Nasal Gouge v-shaped
- 1 Hartmann Nasal Forceps
- 1 S.M.R Hammer
- 1 Frazier Suction Nozel
- 1 Nasal Packing Forceps
- 1 Ballenger Swivel Knife Cvd 4mm

Handwritten signature and initials in blue ink.

- 1 Adson Micro Tissue Forceps
- 1 Luc Turbinate Forceps Medium
- 2 Killin Nasal Speculum
- 1 Heymann Turbinate Scissors Angled
- 3 Thudichum Nasal Speculum
- 1 Jansen Middlition Punch
- 1 Walsham Septum Forceps Straight
- 1 Instruments Box S.S

29-Polypectomy set:



Cold forceps polypectomy (one bite)	Polyp <4 mm in size	Forceps (cup size 2-3 mm)	-
Cold forceps polypectomy (two bite)	Polyp in size of 4-5 mm	Forceps (cup size 2-3 mm) or Jumbo forceps (cup size >3 mm)	-
Hot forceps polypectomy	Polyp <5 mm in size	Hot forceps	Injection needle
Cold snare polypectomy	Polyp in size of 4-5 mm	Mini-oval snare (10-15 mm in diameter)	-
	Polyp in size of 5-10 mm	Mini-oval snare (preferred)	-

Handwritten notes in blue ink:
 Killin
 1. 2. 3.

Hot snare
polypectomy
(small polyp)

Hot snare
polypectomy
(large polyp)

EMR
polypectomy

Polyp >10 mm in size

Polyp in size of 5-10 mm
with flat morphology or
polyps >10 mm in size

Barbed snare (oval snare
with spike or spiral snare) if
the grasp is hard or easily
slipping off

Oval, hexagonal snares

Oval, mini-oval, hexagonal
snares

Detachable snare for
large pedunculated
polyp if needed.

Injection needle

Detachable snare for
large pedunculated
polyp if needed

30-Electro coater Nasal:



Technical specification:

Part of the **Institute of Physics award-winning Ossila Solar Cell Prototyping Platform**, the Ossila Spin Coater offers the ideal spin-coating solution for a busy lab where space is at a premium. **Spin coating** is a widely used and versatile technique for depositing materials onto substrates with accurate and controllable film thicknesses.


The product is compact and highly affordable, allowing everyone to implement spin coating in their research. It offers high-quality coating without the issues of substrate warping. As it does not require a vacuum pump or nitrogen line, it is truly the ultimate plug-and-play spin coater. Its compact size and reduced need for servicing enables you to optimise space in the glovebox or on the bench — without compromising on functionality.

Our Spin Coater is fully compatible with the **Ossila Syringe Pump** and can be purchased as a **package deal**.

[Handwritten signature and scribbles in blue ink]


Designed to reduce substrate warping, improve film quality & eliminate the service requirements necessary with vacuum/nitrogen lines.

31-Nasal bone fracture set:



- Asch septum-straightening forceps
- Walsham septum-straightening forceps
- Boies nasal fracture elevator
- Mayo hemostat with rubber tubing
- Killian nasal septum speculum

32-Tracheostomy set:



Handwritten signature and initials in blue ink.

Simaeco Medical


TRACHEOSTOMY TRAY

Surgical Instruments Sets

- 1-Surgical Knife Handle #3
- 2-Crile Hemostat, 5.5", Cvd.
- 3-Lukens Trachea Retractor, 6.75"
- 4-Trachea Retractor, 6.5", 1 Prong, Sharp
- 5-Waetzlarer Retractor, 5.25", Sharp
- 6-Crile Wood Hoodle Holder, Light
- 7-Skillman Mosquito Forceps, 5", Sur
- 8-Skillman Mosquito Forceps, 5", Cvd
- 9-Allis Forceps, 5 x 6, 6"
- 10-Mayo Poon, 5.5", Cvd.
- 11-Ballenger Sponge Forceps, 7", Cvd
- 12-Backhaus Towel Clamp, 3"
- 13-Jackson Tissue Forceps, 5 x 2, 5"
- 14-Adson Tissue Forceps, 5", Surr. Tips
- 15-Adson Tissue Forceps, 5", Surr. Tips
- 16-Spencer Stitch Scissors, 4.75"
- 17-Mayo Scissors, 6.75", Cvd.
- 18-Motzbaum Scissors, 5.5", Cvd.
- 19-Cushing Vein Retractor, 8"
- 20-Yankauer Suction Tube S.S.




33- Tonsillectomy set:



Handwritten signature and date: 15/11/11


Description	Qty
Yankauer Suct Tube Stainless	1
Mayo Scissors Curved 6 3/4"	1
Rochester-Pean Forceps Cvd 6 1/4"	2
Crile-Wood Needle Holder 7"	1
Allis Tissue Forceps 7 1/2"	1
Backhaus Towel Clamps 5 1/2"	1
White Tonsil Forceps 7 1/2"	2
Schmidt Tonsil Hemostats	2
Jennings Mouth Gag Small	1
Jennings Mouth Gag Medium	1
Jennings Mouth Gag Large	1
Weider Tongue Depressor Small	1
Weider Tongue Depressor Large	1
Love Uvula Retractor	1
Metzenbaum Scissors Cvd 7"	1
Hurd Dissector and Pillar Retractor	1
Eves Tonsil Snare	1
Tyding Tonsil Snare	1
Hartman Tonsil Punch	1
Foerster Sponge Forceps Str Serr 9 1/2"	1

34-Screw Driver:



Handwritten signature and date: 11-1-11

35-Binder of plates:



36-Drill bit:





Technical specification:




Handwritten signature

Handwritten text: 6-1-15


This 9 piece SP lip and spur set from FISCH is precisely ground and sharpened in Austria with the exact geometry needed for fast accurate drilling. Consists of 3,4,5,6,7,8,9 and 10mm bits, supplied in a plastic storage case.

37- Ex-fix apparatus AO set specification:





38-DCP Broad plates:



39-DCP Narrow Plates:

DCP 4.5mm, NARROW



Fixation with 4.5 mm & 6.5 mm screws. All holes are designed as self - compressing holes. All the end holes are designed for cancellous bone screws.


Indication : As a neutralization and tension band plate on the tibia

Profile:


Thickness :	5.0mm
Width :	16.0mm
Hole Space :	16mm & 25mm
Hole :	4 Hole to 18 Hole

Handwritten signature and date: 11/11/11

40-DCP Small:



41-K-wires:



Material: Stainless steel 316L, 316LVM Titanium


Surface Finishing: Powder Coated

Diameter: 1.0mm, 1.2mm, 1.5mm, 1.8mm, 2.0mm, 2.5mm, 3.0mm, 3.5mm


Length: 150mm, 180mm, 200mm, 225mm, 250mm, 300mm

Handwritten signature and initials in blue ink.

42-Wires orthopedic:




43-cortical screw:



[Handwritten signature]

[Handwritten signature]

44-Sponge screw:




Technical specification:

Clearance! Literato - Tool Sponge Brush! Hot Sale! New! Screw Extractor Drill Bits Guide Broken Damaged Bolt Easily Remove Stripped or Damaged Screws Remover Tool Set

[Handwritten signature]

45-Malleolar screw:


4.5mm MALLEOLAR SCREW



Profile:

Thread Diameter :	4.5mm
Drill Bit for Threaded Hole :	3.2mm
Drill Bit for Gliding Hole :	3.0mm
Core Diameter :	3.5mm
Hexagonal Socket :	3.5mm
Head Diameter :	8.0mm
Length :	25mm to 80mm


46-Spongeal screw Partial threaded:




Handwritten signature

Handwritten signature

47-Tubular plates:



48-Pelvic Reconstructive plates:



Reconstruction Plate for Mandible
Straight

Screw Ø : 1.5mm, 2.0mm & 2.5mm
Holes : 8, 10, 12, 16, 20, 24, 30




Reconstruction Plate for Mandible
Angled (Left & Right)

Screw Ø : 1.5mm, 2.0mm & 2.5mm
Holes : 16, 20, 24




Reconstruction Plate for Mandible
Double Angled


Screw Ø : 1.5mm, 2.0mm & 2.5mm
Holes : 4 + 16 + 4, 4 + 20 + 4,
4 + 24 + 4, 4 + 28 + 4



49-T -Plates:



50-sign nail (interlocked set):



[Handwritten signature]

[Handwritten signature]




This set is designed for use in treating patients who need but cannot afford an IM nail and interlocking screw system. They are made from stainless steel that satisfies ASTM implant-grade material specifications acceptable. The orthopedic hardware manufacturing operation is being conducted in accordance with national good manufacturing practices and quality assurance standards.


Handwritten signature and initials in blue ink.

51-Kuntscher Nail set:

Figure 1. Interlocking nail instrumentation and implants.





Complete 6- and 8-mm interlocking nail instrumentation and implant sets.



Individual instruments, including a drill jig (A), drill and tap sleeves (B), short femoral and humeral extensions (C), long tibial extensions (D), 6- and 8-mm reamers (E), a nail insertion tool (F), an attachment screw (G), a trocar (H-I), a depth gauge (J), a nail attachment hex driver (J), drill bits (K), and screwdrivers (L).



52-Orthopedic Drill:



Type: Orthopedic Bone Drill

Kind: Orthopedics

Handwritten signature and scribbles in blue ink.


Group: Adult

Product Type: Surgical Power Drill

Materia: Aluminium Alloy, Stainless Steel

Product Application: Orthopedic Surgery

53-Oscillating bone saw:




Specifications:

- High Quality Korean Medical Equipment
- Product : High Quality Stainless Steel Cannulated Medical Drill, Colors Available : Silver, Black - Upto 1000 RPM
- Product : High Quality Stainless Steel Oscillating Saw, Colors Available : Silver, Black - Upto 1000 RPM

[Handwritten signature and scribbles in blue ink]

54-DHS set:



Description

DHS/DCS Instrument Set

Following Instrument will be there in a SET.

- 1- T-HANDLE WITH QUICK COUPLING
- 2- 135 DEGREE DHS ANGLE GUIDE
- 3- 95 DEGREE DCS ANGLE GUIDE
- 4- VARIABLE ANGLE GUIDE
- 5- DIRECT MEASURING DEVICE

[Handwritten signature and scribbles in blue ink]

6- DHS TRIPLE REAMER

7- DCS TRIPLE REAMER

8- DHS/DCS 12.5MM TAP

9- CENTERING SLEEVE FOR TAP

10- GUIDE SHAFT

11- COUPLING SCREW

12- WRENCH FOR DHS/DCS SCREW



13- CENTERING SLEEVE FOR WRENCH

14- IMPACTOR

15- COUPLING SCREW FOR REMOVAL



16- 2.5MMX230MM THREADED PIN

55-Bipolar set:




Release

1. 2. 3. 4. 5. 6. 7. 8. 9. 10. 11. 12. 13. 14. 15. 16. 17. 18. 19. 20. 21. 22. 23. 24. 25. 26. 27. 28. 29. 30. 31. 32. 33. 34. 35. 36. 37. 38. 39. 40. 41. 42. 43. 44. 45. 46. 47. 48. 49. 50. 51. 52. 53. 54. 55. 56. 57. 58. 59. 60. 61. 62. 63. 64. 65. 66. 67. 68. 69. 70. 71. 72. 73. 74. 75. 76. 77. 78. 79. 80. 81. 82. 83. 84. 85. 86. 87. 88. 89. 90. 91. 92. 93. 94. 95. 96. 97. 98. 99. 100.



56-Locking plate set:




Technical specification:

| | |
|-----------------|---------------------------|
| Material | Stainless Steel, Titanium |
| Brand | Bond well |
| Type | 2 Hole Plate |
| Plate | Locking plates |
| Material Grade | ss 316L |
| Number Of Holes | up to 20 |

Handwritten signature and date: 11.2.15

57-Slit lamp:



58-Massager Machine:

HAN RIVER[®]



Handwritten signature and notes in blue ink.

Specification:

Dimensions:15*16.5*12cm

Weight:870g

Caution:

Do not use the device if you are pregnant or suffering from heart problems or any physical injuries, and consult your physician if you are unsure.

Only use the device one hour after meals.

Keep away from children.

Do not use the device on the same area of body for more than 15 minutes at a time.

Do not use the product on your head.

Do not use the device in an environment that exceeds 37 degrees C or 100 degrees F.

If the mesh cover is torn or the attachment cap is severely damaged, please replace immediately by contacting your local authorized dealer to purchase these parts.

If the user experiences dizziness, nausea, chest pain, or any other abnormal symptoms, stop the workout at once. Consult a physician immediately.

Warning:

Mesh cover can only be used with the Wavy or Flat attachment caps. Do not use mesh cover when roller or filling caps are attached.

Do not use the machine without mesh cover on areas of your body where hair can be entangled with the device (i.e. neck or shoulder, etc.)

NOTE: We will ship the plug adapter to fit your country plug seat.

Package Content:

1 * Body Massager

3* Massage Heads


1 * Massage Cover

1 * User Manual


note:NO carton out box, but with foam.Because it is the outbox is so heavy and we help customer to save the freight fee . Hope your understandings

Item Include

1* Professional Body Massager

Handwritten signature and date in blue ink. The signature is on the left, and the date '11/11/15' is on the right.

59- Interferential unit:



Product Features & Specifications:



- **Channel:** Dual, isolated between channels
- **Mode:** 6 operation programs
- **Waveform:** Symmetrical Balanced Sine Wave
- **Output Current:** Adjustable, 0-33 mA across 500 Ohm load on each channel
- **Beat Frequency:** 2-160 Hz
- **Carrier Frequency:** 4000 Hz (with $\pm 1\%$ tolerance) adjustable on Channel 2
- **Power Source:** 4x 1.5V AA Batteries OR AC/DC adapter (Input: 100-240V AC; Output: 6V, 1A DC)
- **Low Battery Detection:** $4V \pm 0.1V$
- **Treatment Timer:** Continuous, 15, 30, 45, or 60 minutes
- **Usage Timer:** Maximum 999 hours, 59 minutes
- **Automatic Power Off:** If the output has continually remained at zero for 5 minutes, the device will then be automatically shut down
- **Output Configuration**
 - *4-Pads Treatment*
 - Conventional Interferential Current Therapy (IFC) with Channel 1 delivers a fixed 4kHz sine wave, while output of Channel 2 can be varied between 4002-4160 Hz.
 - *2-Pads Treatment*
 - Pre-modulated Interferential Current Therapy (IFC) in which both Channel 1 and 2 will simultaneously deliver an amplitude-modulated output with beat frequency varied between 2-160 Hz.
- **Package includes everything needed to start treatment right away:**
 - 1x IF 4K Device
 - 2x Sets of Specialty IF 4K Lead Wires
 - 4x 1.5" x 1.5" Electrodes

Handwritten blue ink notes:
1. 1.5" x 1.5" Electrodes
2. 2x Sets of Specialty IF 4K Lead Wires
3. 1x IF 4K Device


- 4x AA Batteries
- 1x AC Adapter
- Belt Clip
- Instruction Manual

60- Traction belts:

EASY USAGE:



61-TractionTable:



[Handwritten signature]

[Handwritten signature]

- ### TTET-200 Electric Hi-Lo Traction Platform Specifications:


- Warranty:** One-year warranty

- A laryngoscope is an instrument used to examine the larynx. Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube.
- The Macintosh Laryngoscope Blade . . . The Macintosh is the predominate model of the curved blade types. During laryngoscopy, the tip or beak of the blade is compressed into the angle formed by the base of the tongue and the epiglottis, indirectly raising the epiglottis.
- Laryngoscope Set includes:
 1. Standard Medium handle, Length: 15cm, Weight 6.39oz. (Not Included: Accept: 2 Size C Batteries)
 2. Macintosh Blade 1 (Infant) Overall Length 90mm, Weight: 3.84 cz.




3. Macintosh Blade 2 (Child) Overall Length 100mm, Weight: 3.42 oz.
4. Macintosh Blade 3 (Med. Adult) Overall Length 130mm, Weight: 2.89 oz.
5. Macintosh Blade 4 (Lg. Adult) Overall Length 155mm, Weight: 2.53 oz.
6. Protective Box

63- Abdominal Hysterectomy set:



| <u>Comp NSN</u> | <u>Nomenclature</u> | <u>Qty</u> | <u>Manufacturer</u> |
|-----------------|---------------------|------------|---------------------|
|-----------------|---------------------|------------|---------------------|

| | | | |
|------------------|------------------------|---|----------------|
| 6530-00-197-9228 | Wrapper Ster 54 IN SQ | 2 | |
| 6515-00-149-1103 | Snare Tonsil Tyding DE | 1 | SKLAR |
| 6515-00-299-8737 | Holder Suture Needle H | 2 | Columbiasurgic |



[Handwritten signature]

[Handwritten signature]

| | | | |
|------------------|-------------------------|------|-------------------|
| 6515-00-320-4590 | Forceps Towel Backhaus | 2 | SolwayInc |
| 6515-00-324-5205 | Depressor Tongue Blade | 1 | Dittmar |
| 6515-00-326-2900 | Dissector & Retractor | 1 | EmiltenbergInc |
| 6515-00-334-4300 | Forceps Hemo Rochester | 2 | Columbiasurgic |
| 6515-00-338 2900 | Forceps HemoSchmidt D | 2 | Columbiasurgic |
| 6515-00-344-7880 | handle Surgical Knife | 1 | Propper mfgco |
| 6515-00-345-4300 | Knife Tonsil Douglas 8 | 1 | Miltex Instrument |
| 6515-00-348-5300 | Needle Hypo C26A Tonsil | 1 | Bectondicinson |
| 6515-00-348-5310 | Needle Hypo C25A Tonsil | 1 | Bectondicinson |
| 6515-00-365-7100 | Scissors Tonsil Metzen | 1 | AP |
| 6515-00-386-7600 | Cannula Laryngeal Angu | 1 | Hevesy medical |
| 6515-00-827-3104 | Syringe Hypo GP 3ML CA | 1 | Bectondickinso |
| 6515-00-832-1482 | Forceps Tonsil Seizing | 1 | Starz |
| 6530-01-166-9036 | Tray Instrument CRES 1 | 1 | Starz |
| 7210-00-299-9610 | Towel Hand Cotton 36x1 | 1 | |
| 6530-01-244-2776 | Wrapper Sterlization PG | 0.17 | |

[Handwritten signatures and initials in blue ink]

64-Scale for traction:



Technical specification:

This versatile system provides safe, effective, and easy positioning of the shoulder during all types of arthroscopic or open shoulder surgery performed in the lateral decubitus position. The arm is placed in a shoulder traction and rotation (STaR) sleeve and then connected to the 3 Point Shoulder Traction Tower. Three-point traction with a lateral strap permits ideal shoulder positioning for improved access to the glenohumeral joint. During more routine arthroscopic procedures, single-point traction may be selected at any desired angle of abduction by transferring weights to the third traction cable. The 3-Point Shoulder Distraction System attaches easily to standard OR table clark rails and has color-coded cable ends for easy identification when transferring traction weights.

Handwritten signature in blue ink.

Handwritten signature in blue ink.